

# KD Intra-Articular® Gel

Viscosupplementation Solution for Orthopaedics

GERMANY  
**ALBOMED**



**High comfort in  
osteoarthritis therapy**



# KD Intra-Articular® Gel

## Viscosupplementation Solution for Orthopaedics

KD Intra-Articular® Gel is a high-quality viscoelastic gel, developed, manufactured and sterilized with steam according to the high international manufacturing standards. It contains the required hyaluronic acid with the well-targeted molecular weight. Its integrity and quality are maintained by the double sterile barrier system.

Due to the high purity, biological and non-animal origin of the components, our product does not have any known immunogenic or pathogenic potential. KD Intra-Articular® Gel is highly biocompatible and resorbable.

- Improvement of joint function
- Anti-inflammatory effect
- For all synovial joints
- One-Shot with 120 mg/syringe
- Easy handling
- Wide product range
- Well-targeted molecular weight
- Biofermented origin
- Double sterile barrier system
- CE certified

### ALBOMED GmbH

Hildebrandstrasse 11  
 90592 Schwarzenbruck | Germany  
 Phone: +49 (0) 9183 95 69 82-0  
 info@albomed.eu | www.albomed.eu



Specification	1.0 %	1.0 %	1.6 %	1.6 %	2.2 %	2.5 %	2.5 %
Sodium hyaluronate (mg)	10 mg	20 mg	16 mg	32 mg	44 mg	75 mg	120 mg
Molecular weight* (MDa)	1.6	1.6	1.6	1.6	1.7	2.2	2.2
Viscosity* (mPas)	25 000	25 000	50 000	50 000	300 000	500 000	500 000
Osmolality (mOsm/kg)	270–400	270–400	270–400	270–400	270–400	270–400	270–400
Storage	2–25 °C	2–25 °C	2–25 °C	2–25 °C	2–25 °C	2–25 °C	2–25 °C
pH	6.8–7.4	6.8–7.4	6.8–7.4	6.8–7.4	6.8–7.4	6.8–7.4	6.8–7.4
Volume (ml)	1.0	2.0	1.0	2.0	2.0	3.0	4.8
Shelf life (months)	42	42	42	42	42	42	42

\*mean value after steam sterilization